

Modular Hip Orthosis

The Ortho Europe modular hip orthosis, is constructed using lightweight, high strength materials and provides:

For

- Hip joint stabilisation
- Abduction / adduction control
- Range of motion control



Single hex key to adjust waist size, flexion/extension and adduction /abduction



Anatomically shaped thigh component



Single screw adjustment of abduction/adduction in 7,5 degree increments



Vertical adjustment between hip & thigh components. Interchangeable distal hinge section to accommodate tall patients. Plastic components perforated for ventilation



Multifunctional hinge for simple reliable control of range of motion and abduction / adduction. Simple, precise control of hip flexion / extension in 15 degree increments

Other functions:

- Simple reliable, graduated waist size adjustment. Rotational and lateral adjustment
- Quick release strap with safety lock

Designed for maximum comfort with machine washable padding

Modular Hip Orthosis

Ordering Guidelines



Hip / Waist waits component				
Size	Part Number		Waist Circumference mm (A)	Hip Circumference mm (B)
	Right	Left		
Small	MHH - S - R	MHH - S - L	600 - 820	680 - 900
Medium	MHH - M - R	MHH - M - L	800 - 960	900 - 1100
Large	MHH - L - R	MHH - L - L	940 - 1150	1060 - 1260

Thigh component				
Size	Part Number		Proximal Circumference mm (C)	Distal Circumference mm (D)
	Right	Left		
Small	MHT - S - R	MHT - S - R	440 - 530	340 - 410
Medium	MHT - M - R	MHT - M - R	480 - 600	360 - 450
Large	MHT - L - R	MHT - L - R	570 - 680	420 - 490

"Kinesio" hinge		
Size	Part Number	
Regular	KH - R	Patient < 1.8m
Long	KH - L	Patient > 1.81m